

# Diagnosis of prosthetic mitral valve endocarditis due to *Tropheryma whipplei* Using Next-Generation Sequencing of Plasma



Contact: brian.buggy@aurora.org  
IDWeek Poster 710

Christiaan R. De Vries<sup>1</sup>, Ann T. MacIntyre<sup>1</sup>, Brian P. Buggy<sup>2</sup>

<sup>1</sup> Karius, Inc., Redwood City, CA; <sup>2</sup> Aurora Health Care, Milwaukee, WI

## Background

*Tropheryma whipplei* is a gram-positive bacillus that causes Whipple's disease, a protean multisystem syndrome classically characterized by arthralgias, chronic diarrhea, malabsorption and weight loss. *T. whipplei* infection has a wide spectrum of clinical manifestations including gastrointestinal, musculoskeletal, neurologic, dermatologic and cardiac. Endocarditis has been diagnosed in a small number of patients. Diagnosis is usually accomplished by histopathology on resected valvular tissue or GI tract biopsy with a compatible cardiac course. Next-generation sequencing (NGS) of microbial cell-free DNA (mcfDNA) in plasma offers a rapid, non-invasive means of diagnosis of this rare cause of culture-negative endocarditis.

## Methods

Microbial cell-free DNA was extracted from plasma and NGS was performed by Karius in its CLIA certified/CAP accredited laboratory (Redwood City, CA). Human sequences were removed and the remaining sequences aligned to a curated database of over 1,000 pathogens. Organisms present above a predefined statistical significance threshold are reported and quantified in DNA molecules per microliter (MPM).

## Results

An adult male with prior tissue aortic and mitral valve replacements presented with worsening congestive heart failure. Echocardiography revealed a thickened mitral valve with a small mobile vegetation and severe mitral stenosis. An exhaustive infectious blood culture and serologic evaluation was negative. Karius testing detected *T. whipplei* at 766 MPM within two days of sample receipt (8 days after sample acquisition). The normal range for *T. whipplei* is 0 MPM based on a cohort of 684 healthy individuals. Blood PCR for *T. whipplei* was confirmatory (ARUP Laboratories, Salt Lake City, UT with a turnaround time of 21 days).

## Clinical Parameters of Case of *T. whipplei* infection diagnosed by NGS of plasma mcfDNA

Parameters	Values
Age and Gender	Adult male
Presenting symptoms	Exertional dyspnea
Antecedent symptoms	None
Tmax/Fever at presentation	99.8°F Tmax, otherwise afebrile
Hgb/Hct	8.4/27.5
WBC with %N	8.7 with 83%
Platelets	188,000
PT/PTT	PT 33.2 INR 3.3 No PTT
ESR mm per hr/CRP md per dL	ESR 49 CRP 7.9 mg/dl
Albumin	3.8
Blood culture result	9 sets negative for bacteria
Sites/organ systems involved: Heart	Heart
Imaging results	Pulmonary edema CT chest/abd/pelvis otherwise negative
Abx pretreatment duration prior to Karius Test	Vancomycin/ceftriaxone for 4 days
Choice of antibiotics after Karius Test	Ceftriaxone/moxifloxacin
Karius Test result	<i>Tropheryma whipplei</i> 766 MPM RR 0 MPM (turnaround ~ 2 days)
Other infectious disease testing, result and turnaround time:	<i>T. whipplei</i> to ARUP on 2 occasions (turnaround ~21 days) Negative tests: <i>Histoplasma</i> and <i>Blastomyces</i> antigens, CF and ID antibodies, Fungitell assay, <i>Coxiella</i> serology, <i>Bartonella quintana</i> PCR, <i>Brucella</i> antibodies, Rickettsial antibodies, Blood PCR for CMV, EBV and BKV, <i>Legionella</i> antibody

MPM – molecules of microbial cell-free DNA/microliter; RR – reference range based on the 97.5<sup>th</sup> %ile in a cohort of 684 healthy individuals

## Conclusion

NGS for mcfDNA in plasma offers a rapid, non-invasive method for identifying *T. whipplei*. To our knowledge, this is the first report using NGS of plasma mcfDNA to diagnose prosthetic valve endocarditis due to *T. whipplei*.