

Plazomicin Susceptibility Testing using ETEST[®] MIC for Enterobacterales

Ying Y,¹ Armstrong T,² Blanchard L,³ Denys G,⁴ Halimi D,⁵ Kresken M,⁶ Pompilio M,⁵ and Zambardi G⁵

¹Quest Diagnostics, Lewisville, TX, USA; ²bioMérieux Global Clinical Affairs, St Louis, MO, USA; ³bioMérieux Global Clinical Affairs, Marcy-l'Étoile, France; ⁴Indiana University School of Medicine, Indianapolis, IN, USA; ⁵bioMérieux, La Balme-les-Grottes, France; ⁶Antiinfectives Intelligence GmbH, Reinbach, Germany.

Correspondence: Yun (Claire) Ying
Yun.X.Ying@QuestDiagnostics.com

Introduction

- Plazomicin (PLZ, trade name ZEMDRI[®]) is an aminoglycoside antibiotic that has been shown to be active against *Escherichia coli*, *Klebsiella pneumoniae*, *Proteus mirabilis*, and *Enterobacter cloacae*.
 - PLZ is indicated to treat complicated urinary tract infections (cUTI) in adults who have limited or no alternative treatment options.
- ETEST[®] (bioMérieux SA, Marcy-l'Étoile, France) is a manual, quantitative technique for determining antimicrobial susceptibility of nonfastidious Gram-negative and Gram-positive aerobic bacteria and fastidious bacteria (Figure).
 - The system comprises a predefined antibiotic gradient, which is used to determine the minimum inhibitory concentration (MIC, in µg/mL) of different antimicrobial agents against microorganisms tested on agar media after overnight incubation.¹




Figure. ETEST PLZ test for measuring PLZ MIC, here, on E Coli ATCC[®] 25922[™].

- Here, we report the results of a clinical trial, which compared the performance of ETEST PLZ to the broth microdilution (BMD) reference method.

Methods

- The clinical trial was comprised 4 performance components, which were conducted at the following sites: bioMérieux SA Clinical Affairs Laboratory, Marcy-l'Étoile, France; Indiana University School of Medicine, Indianapolis, IN, USA; Quest Diagnostics at Medfusion, Lewisville, TX, USA; and Antiinfectives Intelligence GmbH, Cologne, Germany.
 - (1 and 2) Clinical and challenge (comparative) studies:** The challenge component included 80 isolates; the clinical component included 518. The combined sets included 598 isolates of species listed in the US Food and Drug Administration (FDA)-approved drug instructions for use (IFU).²
 - (3) Reproducibility study:** 10 isolates were tested in triplicate for 3 days at 3 clinical trial sites.
 - (4) Quality control (QC):** QC was conducted following CLSI QC guidelines³ using the following organisms: *Escherichia coli* ATCC 25922[™], *Enterococcus faecalis* ATCC 29212[™] on ETEST PLZ and BMD; *Pseudomonas aeruginosa* ATCC 27853[™], and *Staphylococcus aureus* ATCC 29213[™] on BMD only.
- Isolates were subcultured on tryptic soy or Columbia agar plates supplemented with 5% sheep blood before testing. After incubation, suspensions of the isolates were prepared in 0.85% saline solution. These suspensions were used to inoculate the BMD or Mueller Hinton agar for the ETEST. Results were read after 16 to 20 hours incubation at 35 °C ± 2 °C in ambient air.
- Results were analyzed using US FDA breakpoints for PLZ: susceptible, ≤2 µg/mL; intermediate, 4 µg/mL; resistant ≥8 µg/mL.⁴ There are no CLSI breakpoints to date.
- Performance was evaluated using US FDA performance criteria: essential agreement (EA) and category agreement (CA) ≥90%, major error rate ≤3.0% and very major error rate ≤2.0%.⁵

Methods (continued)

- EA is defined as** a MIC result obtained with the antimicrobial susceptibility testing system (ASTS) within 1 doubling dilution step; **CA**, agreement of susceptible (S), intermediate (I), or resistant (R) between a MIC test and the reference method; **minor error**, difference in results between the ASTS and the reference method, one result being I and the other is susceptible or resistant; **major error**, reference method result is S and the ASTS result is R; **very major error**, reference method result is R and the ASTS result is S.⁶

Results

- A total of 598 isolates were tested in the comparative studies. The overall performance and MIC distribution are listed in **Tables 1 and 2**. Organisms that tested are listed in **Table 3**.
- Most clinical isolates were susceptible at 0.5 µg/mL PLZ (**Table 2**).

Table 1. Overall Performance of ETEST PLZ for Enterobacterales

EA, % (n)	CA, % (n)	Very major error rate, % (n)	Major error rate, % (n)	Minor error rate, % (n)
99.0 (592/598)	92.8 (555/598)	1.9 (1/53)	0.0 (0/478)	7.0 (42/598)

CA, categorical agreement; EA, essential agreement.

Table 2. Distribution of Test Performance and MICs of Isolates Tested With Each Testing System

ETEST PLZ MIC, µg/mL, and interpretation	Reference method MIC, µg/mL, and interpretation													Total, <i>mode</i>		
	≤0.016 S	0.032 S	0.064 S	0.125 S	0.25 S	0.5 S	1 S	2 S	4 I	8 R	16 R	32 R	64 R		128 R	≥256 R
<0.016 S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0
0.016 S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0
0.032 S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0
0.064 S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0
0.125 S	-	-	-	2	2	-	-	-	-	-	-	-	-	-	-	5
0.25 S	-	-	3	44	10	2	-	-	-	-	-	-	-	-	-	59
0.5 S	-	-	-	75	136	19	2	-	-	-	-	-	-	-	-	232
1 S	-	-	-	1	26	54	12	-	-	-	-	-	-	-	-	93
2 S	-	-	-	-	-	23	52	18	1	-	-	-	-	-	-	94
4 I	-	-	-	-	-	-	15	45	5	-	-	-	-	-	-	65
8 R	-	-	-	-	-	-	-	4	5	2	-	-	-	-	-	11
16 R	-	-	-	-	-	-	-	-	4	3	-	-	-	-	-	7
32 R	-	-	-	-	-	-	-	-	-	2	1	-	-	-	-	3
64 R	-	-	-	-	-	-	-	-	-	-	-	1	-	-	-	0
128 R	-	-	-	-	-	-	-	-	-	-	-	-	1	-	-	0
≥256 R	-	-	-	-	-	-	-	-	-	-	-	-	-	1	29	30
Evaluable	-	-	-	5	122	172	98	81	67	15	7	1	-	-	-	-
Total, <i>mode</i>	0	0	0	5	122	172	98	81	67	15	7	1	0	1	29	-

CA, categorical agreement; EA, essential agreement; I, intermediate; MIC, minimum inhibitory concentration; S, susceptible; R, resistant. Green indicates: EA, CA, and no error; yellow: CA, essential disagreement, and no error; blue: EA, categorical disagreement, and minor error; red: categorical and essential disagreement, and very major error.

Results (continued)

- ETEST PLZ performance was 100% in reproducibility (270/270) and QC studies (81/81 for each QC organism).

Table 3. Distribution of EA, CA, and Errors Among Organisms

Organism	Study Component	EA (of total), n (%)	EA (of evaluable), n (%)	CA, n (%)	Resistant isolates, n	Very major error, n	Major error, n	Minor error, n
<i>C freundii</i>	Challenge	6/6 (100)	6/6 (100)	6/6 (100)	0	-	0	0
	Clinical	53/53 (100)	52/52 (100)	53/53 (100)	1	0	0	0
	Species total	59/59 (100)	58/58 (100)	59/59 (100)	1	0	0	0
<i>C koseri</i>	Challenge	3/3 (100)	3/3 (100)	3/3 (100)	0	-	0	0
	Clinical	30/31 (97)	30/31 (97)	31/31 (100)	0	-	0	0
	Species total	33/34 (97)	33/34 (97)	34/34 (100)	0	-	0	0
<i>E cloacae</i>	Challenge	10/10 (100)	10/10 (100)	10/10 (100)	0	-	0	0
	Clinical	50/50 (100)	48/48 (100)	50/50 (100)	2	0	0	0
	Species total	60/60 (100)	58/58 (100)	60/60 (100)	2	0	0	0
<i>E coli</i>	Challenge	13/13 (100)	8/8 (100)	10/13 (77)	5	0	0	3
	Clinical	63/65 (97)	62/64 (97)	65/65 (100)	1	0	0	0
	Species total	76/78 (97)	70/72 (97)	75/78 (96)	6	0	0	3
<i>K aerogenes</i>	Challenge	6/6 (100)	6/6 (100)	6/6 (100)	1	0	0	0
	Clinical	32/32 (100)	31/31 (100)	32/32 (100)	1	0	0	0
	Species total	38/38 (100)	37/37 (100)	38/38 (100)	2	0	0	0
<i>K oxytoca</i>	Challenge	7/7 (100)	7/7 (100)	6/7 (86)	0	-	0	1
	Clinical	32/32 (100)	31/31 (100)	31/32 (97)	1	0	0	1
	Species total	39/39 (100)	38/38 (100)	37/39 (95)	1	0	0	2
<i>K pneumoniae</i>	Challenge	17/17 (100)	13/13 (100)	16/17 (94)	7	0	0	1
	Clinical	72/72 (100)	63/63 (100)	72/72 (100)	11	0	0	0
	Species total	89/89 (100)	76/76 (100)	88/89 (99)	18	0	0	1
<i>M morgani</i>	Challenge	3/3 (100)	3/3 (100)	2/3 (67)	2	0	0	1
	Clinical	27/28 (96)	27/28 (96)	19/28 (68)	3	0	0	9
	Species total	30/31 (97)	30/31 (97)	21/31 (68)	5	0	0	10
<i>P mirabilis</i>	Challenge	4/4 (100)	4/4 (100)	4/4 (100)	2	0	0	0
	Clinical	59/59 (100)	57/57 (100)	50/59 (85)	6	0	0	9
	Species total	63/63 (100)	61/61 (100)	54/63 (86)	8	0	0	9
<i>P stuartii</i>	Challenge	2/2 (100)	1/1 (100)	2/2 (100)	1	0	0	0
	Clinical	33/33 (100)	30/30 (100)	24/33 (73)	7	0	0	9
	Species total	35/35 (100)	31/31 (100)	26/35 (74)	8	0	0	9
<i>P vulgaris</i>	Challenge	5/5 (100)	5/5 (100)	3/5 (60)	0	-	0	2
	Clinical	29/29 (100)	29/29 (100)	26/29 (90)	1	0	0	3
	Species total	34/34 (100)	34/34 (100)	29/34 (85)	1	0	0	5
<i>S marcescens</i>	Challenge	4/4 (100)	4/4 (100)	4/4 (100)	0	-	0	0
	Clinical	32/34 (94)	32/34 (94)	30/34 (88)	1	1	0	3
	Species total	36/38 (95)	36/38 (95)	34/38 (90)	1	1	0	3
Total	Challenge	80/80 (100)	70/70 (100)	72/80 (90)	18	0	0	8
	Clinical	512/518 (99)	492/498 (99)	483/518 (93)	35	1	0	34
Grand total		592/598 (99)	562/568 (99)	555/598 (93)	53	1	0	42

CA, categorical agreement; EA, essential agreement.

A dash (-) indicates very major error unable to be calculated owing to lack of resistant isolates.

Conclusions

- The clinical performance of the ETEST PLZ test met the US FDA acceptance criteria and was found useful for determining the MIC of Plazomicin for Enterobacterales, including those that produce extended spectrum beta-lactamase (ESBL) or high levels of AmpC β-lactamases; carbapenem-resistant Enterobacteriaceae (CRE, including those that produce metallo-β-lactamases [MBL], *K pneumoniae* carbapenemase [KPC], or oxacillinase-48 [OXA-48]); and aminoglycoside-resistant strains.
- Among the 598 isolates tested, the overall susceptibility to Plazomicin was 80%.
 - The mode MIC was low, at 0.5 µg/mL (susceptible).
- Plazomicin, as an aminoglycoside antibiotic, provides a treatment option for cUTI with multidrug-resistant Enterobacterales.
- The ETEST PLZ test is easy to use and has a quick turnaround time.

References

- ETEST Plazomicin (PLZ). Package insert (IUO). bioMérieux Inc; 2019.
- Drug trials snapshot: ZEMDRI. US Food and Drug Administration. Updated July 12, 2018. Accessed October 8, 2020. <https://www.fda.gov/drugs/drug-approvals-and-databases/drug-trials-snapshot-zemdri>
- M100: Performance Standards for Antimicrobial Susceptibility Testing. 29th ed. Clinical and Laboratory Standards Institute; 2019.
- FDA-recognized antimicrobial susceptibility test interpretive criteria. US Food and Drug Administration. Updated August 27, 2019. Accessed October 8, 2020. <https://www.fda.gov/drugs/development-resources/fda-recognized-antimicrobial-susceptibility-test-interpretive-criteria>
- Antimicrobial susceptibility test (AST) systems - class II special controls guidance for industry and FDA. US Food and Drug Administration. Updated August 28, 2009. Accessed October 8, 2020. <https://www.fda.gov/medical-devices/guidance-documents-medical-devices-and-radiation-emitting-products/antimicrobial-susceptibility-test-ast-systems-class-ii-special-controls-guidance-industry-and-fda>
- M52: Verification of Commercial Microbial Identification and Antimicrobial Susceptibility Testing System. 1st ed. Clinical and Laboratory Standards Institute; 2015.

Acknowledgments

The authors thank the clinical trial staff at the bioMérieux CA Clinical Affairs Laboratory (Marcy CA), Indiana University School of Medicine, Quest Diagnostics at Medfusion (Anita Maharjan, Sarah Fremgen, Tong Zhang, Linda Pham and Lisa Dunnery), and Antiinfectives Intelligence GmbH. We also thank Allison TerBush (Quest, Global Scientific Publications) for assistance preparing the poster.

Disclosure

This work was supported by bioMérieux Clinical Affairs.