

Aram Kadoom¹, Daniel Lockhart¹, Yassar Kadoom¹, Sonia La Fauci¹, James Turner¹, Helen Bennett¹, Paul Jay¹, Sumi Thaker¹, William Mullen¹
¹Momentum Bioscience Ltd, Oxford, UK

Introduction

- Diagnosis of bloodstream infections (BSI) and subsequent treatment with appropriate antimicrobials is dependent upon fast and accurate information about the causative microorganism(s).
- The time taken for blood culture, microbial identification and antimicrobial susceptibility testing (AST), can lead to poor antimicrobial stewardship and patient care.
- Many antimicrobial change decisions are based on the results of a Gram stain, with this being the first result available.
- A rapid test, which can confirm BSI and characterize the causative pathogen(s), would improve antimicrobial stewardship and patient care.

Method

- SepsisSTAT[®] is a rapid direct-from-blood molecular test, developed by Momentum Bioscience Ltd, for the detection of BSI, with a time-to-result of < 4 hours.
- Uses whole blood to detect viable microorganisms whilst also providing molecular characterization (Gram status and genus of key infectious organisms).
- Microorganisms are extracted from the sample through capture on magnetic microbeads, followed by Enzymatic Template Generation and Amplification (ETGA[®]) for ultra-sensitive, universal detection of viable bacterial and fungal species.
- Simultaneously, molecular characterization (Molecular ID) also provides genus/species identification based on total microbial DNA present (from viable cells and cell-free DNA).
- The detection limits of SepsisSTAT[®] were evaluated for a broad panel of microorganisms, representing 80.4% of BSI reported to Public Health England (2018 report).

Results




Figure 1: ETGA results for 16 microorganism panel (n=5). The boxes show the upper/lower quartiles, with the line representing the median cfu/mL. The lines (whiskers) show the lowest and highest detected cfu/mL, with dots representing outliers.




Figure 2: Molecular ID results for 16 microorganism panel (n=5). The boxes show the upper/lower quartiles, with the line representing the median cfu/mL. The lines (whiskers) show the lowest and highest detected cfu/mL, with dots representing outliers. Note that *S. aureus* has two SepsisID detection channels: *Staphylococci* and *Coagulase*

Discussion

- This study was carried out with five biological replicates of each of the microorganisms.
- Of the 16 microorganisms tested, all but two had an ETGA median limit of detection (LoD) of < 10 cfu/mL, and of these, five organisms were detected at < 1 cfu/mL. These results demonstrate that 77.4% of monomicrobial infections are detected at < 10 cfu/mL.
- For Molecular ID, 13 of 16 had a median LoD of < 100 cfu/mL, with six organisms being detected at < 10 cfu/mL. This shows that 76.9% of monomicrobial infections are detected at < 100 cfu/mL. Importantly, total microbial DNA in BSI has been shown to have a typical range of 100-10,000 cfu equivalents/mL (Bacconi *et al.*).
- For ETGA, of the five ESKAPE organisms tested (*A. baumannii* not tested), four of these are detected at levels < 10 cfu/mL, with only *E. faecium* detected at 260 cfu/mL.
- For Molecular ID, four organisms are also detected at levels < 100 cfu/mL, with only *P. aeruginosa* detected at 336 cfu/mL.

References: Bacconi, Andrea, *et al.* (2014). "Improved sensitivity for molecular detection of bacterial and Candida infections in blood." *Journal of clinical microbiology* 52.9: 3164-3174.

Table: Median detections for ETGA and Molecular ID for the 16 microorganisms tested. Key shows cfu/mL blood detection levels and number of microorganisms detected at those levels

Organism	PHE Rank	% Infections	Cumulative Infection	ETGA	Molecular ID
<i>S. epidermidis</i>	1	25.7	25.7	8.4	16.8
<i>E. coli</i>	2	24.9	50.6	5.8	64.2
<i>S. aureus</i> (*Coagulase)	3	7.7	58.3	0.1	3.0/ *3.6
<i>K. pneumoniae</i>	4	4	62.3	3.5	34.8
<i>S. pneumoniae</i>	5	3.7	66	4.6	47.8
<i>P. aeruginosa</i>	6	2.2	68.2	0.5	336.0
<i>S. pyogenes</i>	7	2	70.2	2.4	36.0
<i>E. faecalis</i>	8	1.8	72	0.1	3.6
<i>P. mirabilis</i>	9	1.8	73.8	0.5	5.0
<i>S. agalactiae</i>	10	1.6	75.4	50.8	59.0
<i>E. faecium</i>	11	1.4	76.8	260.0	30.2
<i>S. dysgalactiae</i>	12	1	77.8	2.9	2.7
<i>M. luteus</i>	14	0.9	78.7	8.3	312.0
<i>E. cloacae</i>	15	0.8	79.5	1.0	3.4
<i>S. marcescens</i>	18	0.5	80	0.6	62.6
<i>C. albicans</i>	25	0.4	80.4	6.6	368.0

Key:

	ETGA	Molecular ID
<1000	16	16
<100	15	13
<10	14	6
<1	5	0

Conclusion

- SepsisSTAT[®] results demonstrate the sensitivity of universal detection of viable microorganisms.
- SepsisSTAT[®] detects clinically relevant microorganisms at low levels with sensitive and specific detection.
- Current development aims to shorten the time-to-result to < 3 hours.
- Studies in a clinical setting with seek to further demonstrate the efficacy and potential to impact patient care.